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
A History of Automated Ventilation

Despite years of research in mechanical ventilation many settings remained to be set by hand.

Here is when the ventilators integrated closed-loop systems.

1992

PAV+

monitors flow and volume
inspiratory assist in proportion to patient's effort

1998

ASV

automated selection of V_T and RR
according to the least WOB (Otis)

2002

NAVA

monitors diaphragm activity
inspiratory assist in proportion to diaphragm activity

2006

SmartCare

automated weaning mode
monitors V_T , RR, and $etCO_2$
reduction of PS and performs SBT

2011

INTELLIVENT-ASV

automated selection of V_T and RR
according to the least WOB (Otis) and FOB (Mead)
monitors V_T , RR, $etCO_2$ and SpO_2
automated titrations of AMV, PEEP and FiO_2
reduction of AMV via reduction of PS, performs SBT

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References

1. MacIntyre N, Rackley C, Khusid F (2021) Fifty years of mechanical ventilation-1970s to 2020. *Crit Care Med* 49(4):558–574
2. Damiani E, Adrario E, Girardis M, Romano R, Pelaia P, Singer M et al (2014) Arterial hyperoxia and mortality in critically ill patients: a systematic review and meta-analysis. *Crit Care* 18(6):711
3. Amato MB, Meade MO, Slutsky AS, Brochard L, Costa EL, Schoenfeld DA et al (2015) Driving pressure and survival in the acute respiratory distress syndrome. *N Engl J Med* 372(8):747–755CAS
4. Urner M, Jüni P, Hansen B, Wettstein MS, Ferguson ND, Fan E (2020) Time-varying intensity of mechanical ventilation and mortality in patients with acute respiratory failure: a registry-based, prospective cohort study. *Lancet Respir Med* 8(9):905–913
5. Angus DC, Kelley MA, Schmitz RJ, White A, Popovich J Jr (2000) Caring for the critically ill patient. Current and projected workforce requirements for care of the critically ill and patients with pulmonary disease: can we meet the requirements of an aging population? *JAMA* 284(21):2762–2770CAS
6. British Medical Association (2021) Medical staffing in England: a defining moment for doctors and patients. BMA House, London
7. Mamdani M, Slutsky AS (2021) Artificial intelligence in intensive care medicine. *Intensive Care Med* 47(2):147–149
8. Botta M, Wenstedt EFE, Tsonas AM, Buiteman-Kruizinga LA, van Meenen DMP, Korsten HHM et al (2021) Effectiveness, safety and efficacy of INTELLiVENT-adaptive support ventilation, a closed-loop ventilation mode for use in ICU patients—a systematic review. *Expert Rev Respir Med* 15(11):1403–1413CAS
9. De Bie AJR, Neto AS, van Meenen DM, Bouwman AR, Roos AN, Lameijer JR et al (2020) Fully automated postoperative ventilation in cardiac surgery patients: a randomised clinical trial. *Br J Anaesth* 125(5):739–749
10. Buiteman-Kruizinga LA, Mkadmi HE, Schultz MJ, Tangkau PL, van der Heiden PLJ (2021) Comparison of mechanical power during adaptive support ventilation versus nonautomated pressure-controlled ventilation-a pilot study. *Crit Care Explor* 3(2):e0335
11. Buiteman-Kruizinga LA, Mkadmi HE, Serpa Neto A, Kruizinga MD, Botta M, Schultz MJ et al (2021) Effect of INTELLiVENT-ASV versus conventional ventilation on ventilation intensity in patients with COVID-19 ARDS-An observational study. *J Clin Med* 10(22):5409CAS

12. Georgopoulos D, Xirouchaki N, Tzanakis N, Younes M (2016) Data on respiratory variables in critically ill patients with acute respiratory failure placed on proportional assist ventilation with load adjustable gain factors (PAV+). *Data Brief* 8:484–493
13. Burns KE, Lellouche F, Nisenbaum R, Lessard MR, Friedrich JO (2014) Automated weaning and SBT systems versus non-automated weaning strategies for weaning time in invasively ventilated critically ill adults. *Cochrane Database Syst Rev* 2014(9):Cd008638 Central
14. Kampolis CF, Mermiri M, Mavrovounis G, Koutsoukou A, Loukeri AA, Pantazopoulos I (2021) Comparison of advanced closed-loop ventilation modes with pressure support ventilation for weaning from mechanical ventilation in adults: a systematic review and meta-analysis. *J Crit Care* 68:1–9
15. Bialais E, Wittebole X, Vignaux L, Roeseler J, Wysocki M, Meyer J et al (2016) Closed-loop ventilation mode (IntelliVent®-ASV) in intensive care unit: a randomized trial. *Minerva Anesthesiol* 82(6):657–668