Edition: BP 2025 (Ph. Eur. 11.6 update)

Kanamycin Sulfate

General Notices

Kanamycin Sulphate

(Kanamycin Monosulfate, Ph. Eur. monograph 0032)

C₁₈H₃₈N₄O₁₅S,H₂O 601 5965-95-7

Action and use

Aminoglycoside antibacterial.

Ph Eur

DEFINITION

 $6-O-(3-Amino-3-deoxy-\alpha-D-glucopyranosyl)-4-O-(6-amino-6-deoxy-\alpha-D-glucopyranosyl)-2-deoxy-D-streptamine sulfate monohydrate.$

Antimicrobial substance produced by the growth of certain strains of *Streptomyces kanamyceticus*.

Content

Minimum 750 IU/mg (dried substance).

PRODUCTION

It is produced by methods of manufacture designed to eliminate or minimise substances lowering blood pressure.

CHARACTERS

Appearance

White or almost white, crystalline powder.

Solubility

Freely soluble in water, practically insoluble in acetone and in ethanol (96 per cent).

IDENTIFICATION

A. Thin-layer chromatography (2.2.27).

Test solution Dissolve 10 mg of the substance to be examined in <u>water R</u> and dilute to 10 mL with the same solvent.

Reference solution (a) Dissolve 10 mg of <u>kanamycin monosulfate CRS</u> in <u>water R</u> and dilute to 10 mL with the same solvent.

Reference solution (b) Dissolve 10 mg of <u>kanamycin monosulfate CRS</u>, 10 mg of <u>neomycin sulfate CRS</u> and 10 mg of <u>streptomycin sulfate for identification CRS</u> in <u>water R</u> and dilute to 10 mL with the same solvent.

Plate Suitable plate coated with a 0.75 mm layer of a mixture prepared as follows: mix 0.3 g of <u>carbomer R</u> with 240 mL of <u>water R</u> and allow to stand, with moderate shaking, for 1 h; adjust to pH 7 by the gradual addition, with continuous shaking, of <u>dilute sodium hydroxide solution R</u> and add 30 g of <u>silica gel H R</u>.

Pretreatment Heat the plate at 110 °C for 1 h, allow to cool and use immediately.

Mobile phase 70 g/L solution of potassium dihydrogen phosphate R.

Application 10 µL.

Development Over a path of 12 cm.

Drying In a current of warm air.

Detection Spray with a mixture of equal volumes of a 2 g/L solution of <u>1,3-dihydroxynaphthalene R</u> in <u>ethanol (96 per cent) R</u> and a 460 g/L solution of <u>sulfuric acid R</u>. Heat at 150 °C for 5 min to 10 min.

System suitability Reference solution (b):

— the chromatogram shows 3 clearly separated spots.

Results The principal spot in the chromatogram obtained with the test solution is similar in position, colour and size to the principal spot in the chromatogram obtained with reference solution (a).

- B. Dissolve 0.5 g in 10 mL of <u>water R</u>. Add 10 mL of <u>picric acid solution R</u>. Initiate crystallisation if necessary by scratching the wall of the tube with a glass rod and allow to stand. Collect the crystals, wash with 20 mL of <u>water R</u> and filter. Dry at 100 °C. The crystals melt (<u>2.2.14</u>) at about 235 °C, with decomposition.
- C. Dissolve about 50 mg in 2 mL of <u>water R</u>. Add 1 mL of a 10 g/L solution of <u>ninhydrin R</u> and heat for a few minutes on a water-bath. A violet colour develops.
- D. It gives the reactions of sulfates (2.3.1).

https://nhathuocngocanh.com/bp

TESTS

Solution S

Dissolve 0.20 g in *carbon dioxide-free water R* and dilute to 20.0 mL with the same solvent.

pH (2.2.3)

6.5 to 8.5 for solution S.

Specific optical rotation (2.2.7)

+ 112 to + 123 (dried substance), determined on solution S.

Kanamycin B

Thin-layer chromatography (2.2.27) as described under Identification A with the following modifications.

Test solution Dissolve 0.1 g of the substance to be examined in <u>water R</u> and dilute to 20 mL with the same solvent.

Reference solution Dissolve 4 mg of <u>kanamycin B sulfate CRS</u> in <u>water R</u> and dilute to 20 mL with the same solvent.

Application 4 µL.

Detection Spray with <u>ninhydrin and stannous chloride reagent R</u>. Heat at 110 °C for 15 min.

Limit:

— *kanamycin B*: any spot corresponding to kanamycin B in the chromatogram obtained with the test solution is not more intense than the spot in the chromatogram obtained with the reference solution (4.0 per cent).

Loss on drying (2.2.32)

Maximum 1.5 per cent, determined on 1.000 g by drying at 60 °C at a pressure not exceeding 670 Pa for 3 h.

Sulfated ash (2.4.14)

Maximum 0.5 per cent, determined on 1.0 g.

Sulfate

15.0 per cent to 17.0 per cent of sulfate (dried substance).

Dissolve 0.250 g in 100 mL of <u>water R</u> and adjust the solution to pH 11 with <u>concentrated ammonia R</u>. Add 10.0 mL of <u>0.1 M barium chloride</u> and about 0.5 mg of <u>phthalein purple R</u>. Titrate with <u>0.1 M sodium edetate</u> adding 50 mL of <u>ethanol (96 per cent) R</u> when the colour of the solution begins to change and continue the titration until the violet-blue colour disappears.

1 mL of 0.1 M barium chloride is equivalent to 9.606 mg of SO₄.

https://nhathuocngocanh.com/bp

Pyrogens (2.6.8)

If intended for use in the manufacture of parenteral preparations without a further appropriate procedure for the removal of pyrogens, it complies with the test for pyrogens. Inject per kilogram of the rabbit's mass 1 mL of a 10 mg/mL solution of the substance to be examined in *water for injections R*.

ASSAY

Carry out the microbiological assay of antibiotics (<u>2.7.2</u>). Use <u>kanamycin monosulfate CRS</u> as the reference substance.

STORAGE

If the substance is sterile, store in a sterile, tamper-evident container.

Ph Eur