# **Quality standards**

Edition: BP 2025 (Ph. Eur. 11.6 update)

# Sterile Linen Thread in Distributor

**General Notices** 

(Linen Thread, Sterile, in Distributor for Veterinary Use, Ph. Eur. monograph 0608)

Ph Eur

### **DEFINITION**

Sterile linen thread in distributor for veterinary use consists of the pericyclic fibres of the stem of *Linum usitatissimum* L. The elementary fibres, 2.5-5 cm long, are assembled in bundles 30-80 cm long and spun into continuous lengths of suitable diameter. The thread may be creamy-white or may be coloured with colouring matter authorised by the competer authority. The thread is then sterilised.

### **IDENTIFICATION**

- A. Dissect the end of a thread, using a needle or fine tweezers, to isolate a few individual fibres. Examined under a microscope, the fibres are seen to be 12-31 µm wide and, along the greater part of their length, have thick walls, sometimes marked with fine longitudinal striations, and a narrow lumen. The fibres gradually narrow to a long, fine point. Sometimes there are unilateral swellings with transverse lines.
- B. Infrared absorption spectrophotometry (2.2.24). Examine by attenuated total reflection (ATR).

Absorption maxima and intensities  $3326 \pm 10 \text{ cm}^{-1}$  (medium);  $2911 \pm 12 \text{ cm}^{-1}$  (medium);  $1645 \pm 10 \text{ cm}^{-1}$  (weak);  $1426 \pm 5 \text{ cm}^{-1}$  (medium);  $1315 \pm 3 \text{ cm}^{-1}$  (medium);  $1154 \pm 7 \text{ cm}^{-1}$  (medium);  $1104 \pm 5 \text{ cm}^{-1}$  (medium);  $1050 \pm 5 \text{ cm}^{-1}$  (strong);  $1026 \pm 7 \text{ cm}^{-1}$  (strong).

## **TESTS**

It complies with the tests prescribed in the monograph <u>Strands, sterile non-absorbable, in distributor for veterinary use (0605)</u>.

If stored in a dry state, expose to an atmosphere with a relative humidity of  $65 \pm 5$  per cent at  $20 \pm 2$  °C for 4 h immediate before measuring the diameter and for the determination of minimum breaking load immerse in <u>water R</u> at room temperature for 30 min immediately before carrying out the test.

## **STORAGE**

See the monograph Strands, sterile non-absorbable, in distributor for veterinary use (0605).

### **LABELLING**

See the monograph <u>Strands, sterile non-absorbable, in distributor for veterinary use (0605)</u>.

